



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL
OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY

APPLICATION #: OPM-0486-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: ☒ New ☐ Renewal ☐ Update to Pre-CBC 2013 OPA Number: _____

Manufacturer Information

Manufacturer: Belimed

Manufacturer's Technical Representative: Jay Upchurch

Mailing Address: 8351 Palmetto Commerce Parkway, Suite 101, Ladson, SC 29456

Telephone: 843-216-7424; ext. 133 Email: Jay.Upchurch@belimed.com

Product Information

Product Name: MST-V 6 Series VS1 and VS2 Sterilizers

Product Type: Sterilizers

Product Model Number: 6-0-6 VS1, 6-0-6 VS1E, 6-0-6 VS2, 6-0-6 VS2E, 6-0-9 VS1, 6-0-9 VS1E, 6-0-9 VS2, 6-0-9 VS2E, 6-0-12 VS1, 6-0-12 VS1E, 6-0-12VS2, 6-0-12 VS2E

General Description: MST-V 6 series VS1 and VS2 Series of steam sterilizers and designed for efficiency and superior workflow. An advanced vacuum pump and cooling system ensures minimal water consumption.

Applicant Information

Applicant Company Name: ISAT Seismic Bracing

Contact Person: William V Joerger

Mailing Address: 1020 Crews Road, Suite Q, Matthews NC 28105

Telephone: 510-714-0216 Email: wvjoerger@isatsb.com


I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2016.

Signature of Applicant:  Date: April 23, 2018

Title: Principal Structural Engineer Company Name: ISAT Seismic Bracing

"Access to Safe, Quality Healthcare Environments that Meet California's

STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY
OSH-FD-700 (REV 12/16/15)



OSHPD

Page 1 of 2



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations

Company Name: ISAT Seismic Bracing

Name: William V Joerger California License Number: SE 4545

Mailing Address: 1020 Crews Road, Suite Q, Matthews NC 28105

Telephone: 510-714-0216 Email: wvjoerger@isatsb.com

OSHPD Special Seismic Certification Preapproval (OSP)

- ☐ Special Seismic Certification is preapproved under OSP-
(Separate application for OSP is required)
- ☐ Special Seismic Certification is not preapproved

Certification Method(s)

- ☐ Testing in accordance with: ☐ ICC-ES AC156 ☐ FM 1950-16
- ☐ Other* (Please Specify): _____

*Use of criteria other than those adopted by the California Building Standards Code, 2016 (CBSC 2016) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2016 may be used when approved by OSHPD prior to testing.

- ☒ Analysis
- ☐ Experience Data
- ☐ Combination of Testing, Analysis, and/or Experience Data (Please Specify): _____

List of Attachments Supporting the Manufacturer's Certification

- ☐ Test Report ☒ Drawings ☒ Calculations ☐ Manufacturer's Catalog
- ☐ Other(s) (Please Specify): _____

OFFICE USE ONLY – OSHPD APPROVAL VALID FOR CBC 2016 & ALL PRE-2016 CODE BASED PROJECTS


Signature: Sonia Eliseo Date: 8/9/2019

Print Name: Sonia Eliseo

Title: Senior Structural Engineer

Condition of Approval (if applicable): _____

"Access to Safe, Quality Healthcare Environments that Meet California's





TOTAL SUPPORT

Innovation • Engineering • BIM • Fabrication

A Division of Tomarco Contractor Specialties

Submittal Documents

OSHPD OPM-0486-13

DATE: 08/09/2019

ATTACHMENT OPM DRAWINGS MST-V VS-1 AND VS-2 STERILIZERS

BELIMED

ISAT
1020 Crews Road Suite Q
Matthews, N.C. 28105
704-841-4080



WV 19 Jun 19

FILE NO.: CLT-0117-009

"Empowered by Experience"

REV 3

OSHPD OPM-0486-13 DWG - i



OSHPD OPM-0486-13

DRAWING INDEX

DRAWING INDEX

Cover Page

p i

Index Page

p ii

Drawings for OPM-0486-13

General Notes	p 1
Attachment Notes	p 2
Belimed MST-V 6-0-X VS1 and VS-2 Layout Plan	p 3
Equipment Anchorage at Grade	p 4
Equipment Anchorage at Elevated Slabs	p 5
Dimensions and Attachment Forces	p 6
Belimed MST-V 6-0-6 ELD at Grade $S_{DS} = 2.5$	p 7
Belimed MST-V 6-0-9 ELD at Grade $S_{DS} = 2.5$	p 8
Belimed MST-V 6-0-12 ELD at Grade $S_{DS} = 2.5$	p 9
Brackets	p 10
ELD at Grade Dimensions and Attachment Forces	p 11

OSHPD OPM-0486-13

MANUFACTURE: BELIMED

EQUIPMENT TYPE: WASHERS / STERILIZERS

GENERAL NOTES:

1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2016. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2016.
2. SEISMIC CRITERIA USED: $S_{DS} = 2.50$ $I_p = 1.5$ $a_p = 1.0$ $R_p = 1.5$ (OTHER EQUIPMENT). FOR $z/h = 0$ $F_pH = 1.13$ AND FOR $z/h \leq 1.0$ $F_pH = 3.00$ AND $F_pV = 0.50$. FOR "ELD" MODELS AT GRADE LOCATIONS WHERE (4) ANCHORS ARE USED $S_{DS} = 2.25$, $F_pH = 1.01$ AND $F_pV = 0.45$. FOR "ELD" MODELS ON ELEVATED SLABS $S_{DS} = 2.50$, $F_pH = 3.00$ AND $F_pV = 0.50$.
3. SUPPORT AND ATTACHMENT FORCES ARE DETERMINED USING ASCE 7-10 CHAPTER 13 "SEISMIC DESIGN REQUIREMENTS FOR NONSTRUCTURAL COMPONENTS". AN OVERSTRENGTH FACTOR $\Omega_0 = 1.5$ IS USED FOR CONCRETE MATERIALS PER ASCE 7-10 SUPPLEMENT 1 TABLE 13.6-1. LOADS SHOWN ARE STRENGTH DESIGN LOADS PER CBC 2016 SECTION 1605A.2.
4. THIS PREAPPROVAL COVERS ONLY THE SUPPORTS AND ATTACHMENTS OF THE EQUIPMENT TO THE STRUCTURE.
5. STEEL MATERIALS: PLATE 304 STAINLESS STEEL BY MANUFACTIRER ($F_y = 31.2$ KSI), ALL THREAD ROD ASTM F593 CW1 304 STAINLESS STEEL FOR BOLTS AT ELEVATED FLOORS, ANGLE ASTM A36 ($F_u = 58$ KSI) HOT DIPPED GALVANIZED.
6. CONCRETE SLABS:
 - a. FOR ELEVATED SOLID CONCRETE SLABS: 6" THICKNESS OF NORMAL WEIGHT CONCRETE WITH 3000 PSI MINIMUM STRENGTH.
 - b. METAL DECK: 3" DEEP COMPOSITE STEEL DECK, 20 GAGE MINIMUM, 4 1/2 INCH MINIMUM BOTTOM FLUTE WIDTH AND MINIMUM FLUTE SPACING OF 12", WITH 3 1/4 INCH SAND LIGHT WEIGHT CONCRETE CONCRETE COVER AT 3000 PSI MINIMUM
 - c. FOR SLAB ON GRADE: 6" THICKNESS NORMAL WEIGHT CONCRETE AT 4000 PSI MINIMUM STRENGTH.
7. POST-INSTALLED CONCRETE ANCHORS: HILTI KWIK BOLT TZ (ESR-1917) STAINLESS STEEL 5/8" DIAMETER x 4" EFFECTIVE EMBEDMENT. DRILL 11/16" DIAMETER x 4 3/4" INCH DEEP HOLE; CLEAN HOLE PER MANUFACTURER'S INSTRUCTIONS. FOR ANCHORS INTO THE SOFFITT OF THE METAL DECK USE KWIK BOLT TZ (ESR-1917) 1/2" DIAMETER x 4" MIN. HOLE DEPTH FOR 3 1/4" EFFECTIVE EMBEDMENT WITH 1" MAXIMUM OFFSET FROM THE CENTER OF THE METAL DECK FLUTE. BOTH ARE SUPPLIED BY INSTALLATION CONTRACTOR.
8. EXCERCISE DUE CARE WHEN DRILLING POST-INSTALLED ANCHORS TO AVOID DAMAGING CONCRETE REINFORCEMENT OR


BY: Sonia Eliseo

DATE: 08/09/2019

OPM-0486-13 BELIMED MST-V VS1 AND VS2 STERILIZER GENERAL NOTES



**A Division of Tomarco Contractor Specialties
International Seismic Application Technology**
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com



DRAWN BY: WVJ
DATE: 4/23/18

REVISED BY: WVJ
DATE: 07/19/19
REV NO: 3

SCALE
N.T.S.

PAGE
GEN NOTES

OSHPD OPM-0486-13 DWG - 1

OSHPD OPM-0486-13

MANUFACTURE: BELIMED

EQUIPMENT TYPE: WASHERS / STERILIZERS

ATTACHMENT NOTES:

1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE CBC 2016. THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE CBC 2016.
2. PERIODIC SPECIAL INSPECTION PER CBC 2016 SECTION 1705A AND TABLE 1705A.3 INCLUDING VERIFICATION OF ANCHOR TYPE, ANCHOR DIMENSIONS, CONCRETE TYPE, CONCRETE COMPRESSIVE STRENGTH, ANCHOR SPACING, EDGE DISTANCES, CONCRETE MEMBER THICKNESS, HOLE DIMENSIONS, ANCHOR EMBEDMENT AND ADHERENCE TO THE MANUFACTURER'S PRINTED INSTALLATION INSTRUCTIONS. IN ADDITION, FOLLOW THE PROVISIONS OF THE 2016 CALIFORNIA BUILDING CODE SECTION 1910A.5.5. TORQUE CONTROLLED POST-INSTALLED ANCHORS - TEST USING A CALIBRATED TORQUE WRENCH; 60 FOOT-POUNDS TORQUE FOR 5/8" ANCHORS AND 40 FOOT-POUNDS FOR 1/2" ANCHORS SHALL BE ACHIEVED WITHIN ONE-HALF TURN OF THE NUT. TEST 50% OF THE ANCHORS FOR EACH PIECE OF EQUIPMENT. IF ANY ANCHOR FAILS TEST ALL ANCHORS. REPORT OF TEST RESULTS ARE TO BE SUBMITTED TO THE ENFORCEMENT AGENCY. THE SEOR SHALL PROVIDE REMEDIAL ANCHORAGE DETAILS IN THE EVENT THAT AN ANCHOR FAILS TO MEET THE TEST REQUIREMENTS. FOR THROUGH BOLTS MARK THE NUT LOCATION AT SNUG TIGHT CONDITION. INSPECTOR IS TO VERIFY 3/4 TURN.
3. STRENGTH DESIGN WAS USED FOR ANCHOR FORCE CALCULATIONS INCLUDING Ω_0 PER ACI 318-14 WHERE REQUIRED FOR ATTACHMENT TO CONCRETE.
4. PROVIDE FOR FULL THREAD ENGAGEMENT OF THE NUT AND WASHER.


RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD

1. CONFIRM THE MATERIAL PROPERTIES AND THICKNESS OF THE CONCRETE SLAB TO WHICH THE EQUIPMENT IS ATTACHED MEETS THE REQUIREMENTS OF THIS OPM.
2. PROVIDE A PLAN FOR INSPECTION OF SUPPORTS AND ATTACHMENTS AND VERIFY ITS IMPLEMENTATION.
3. CONFIRM THE SPECIFIED MINIMUM CONCRETE EDGE DISTANCES ARE MAINTAINED BASED ON THE ACTUAL EQUIPMENT LOCATION. VERIFY THAT EXISTING OR NEW ANCHORS ARE AN ADEQUATE DISTANCE FROM THIS UNIT'S ATTACHMENT.
4. VERIFY THAT THE EXISTING STRUCTURE IS ADEQUATE FOR THE IMPOSED DEAD, LATERAL AND TENSION FORCES SHOWN IN ADDITION TO ALL OTHER LOADS.
5. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH CBC 2016 AND WITH THE OPM-0486-13 DETAILS INCLUDING MATERIALS AND DIMENSIONS OF THE SUPPORT WHERE THE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN.
6. VERIFY THAT THE PROJECT SPECIFIC S_{DS} AND z/h VALUES RESULT IN SEISMIC FORCES (E_h AND E_v) DO NOT EXCEED THE VALUES SHOWN IN THESE DETAILS.

OPM-0486-13 BELIMED MST-V VS-1 AND VS-2 V446 STERILIZER ATTACHMENT NOTES



**A Division of Tomarco Contractor Specialties
International Seismic Application Technology**
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com




DRAWN BY: WVJ
DATE: 04/23/18

REVISED BY: WVJ
DATE: 07/19/19
REV NO: 3

SCALE
N.T.S.

PAGE
ATTACHMENT


OSHPD OPM-0486-13 DWG - 2



OPM-0486-13 BELIMED MST-V 6-0-X VS1 AND MST-V 6-0-X VS2
 ANCHOR LAYOUT AT GRADE AND ELEVATED SLABS WHERE $S_{DS} = 2.50$
 AND MST-V 6-0-X ELD MODELS AT GRADE WHERE $S_{DS} \leq 2.25$



**A Division of Tomarco Contractor Specialties
 International Seismic Application Technology**
 1020 Crews Road, Suite Q, Matthews, NC 28105
 704-841-4080 www.isatsb.com




DRAWN BY: WVJ
DATE: 04/23/18



REVISED BY: WVJ
DATE: 07/19/19
REV NO: 3

SCALE
 N.T.S.


PAGE
 V60X LAYOUT

OSHPD OPM-0486-13 DWG - 3



 <p>A Division of Tomarco Contractor Specialties International Seismic Application Technology 1020 Crews Road, Suite Q, Matthews, NC 28105 704-841-4080 www.isatsb.com</p>		
	DRAWN BY: WVJ DATE: 04/23/19	
	REVISED BY: WVJ DATE: 07/19/19 REV NO: 3	
	<table border="1"> <tr> <td>SCALE N.T.S.</td> <td>PAGE GRADE</td> </tr> </table>	SCALE N.T.S.
SCALE N.T.S.	PAGE GRADE	

OSHPD OPM-0486-13 DWG - 4




NOTES:

1. DRILL 11/16" HOLE IN THE CONCRETE SLAB TO RECEIVE THE 5/8" ASTM F593 CW1 304 STAINLESS STEEL RODS WITH MATCHING WASHERS AND NUTS.
2. SUPPLEMENTAL ANGLE SHALL BE L3x3x3/8 ASTM A36. DRILL 11/16" HOLES IN ANGLE FOR EQUIPMENT BOLTS. FOR LOCATIONS WHERE A TOP NUT CANNOT BE INSTALLED, DRILL AND TAP THE ANGLE TO RECEIVE THE 5/8" EQUIPMENT ANCHORAGE BOLT.
3. ANCHOR INTO BOTTOM OF SLAB SHALL BE HILTI KWIK BOLT TZ 1/2", SEE NOTE 7 ON PAGE "GEN NOTES".

OPM-0486-13 BELIMED MST-V 6-0-X STERILIZER
ATTACHMENT DETAILS AT ELEVATED SLABS



**A Division of Tomarco Contractor Specialties
International Seismic Application Technology**
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com



DRAWN BY: WVJ	
DATE: 04/23/18	
REVISED BY: WVJ	
DATE: 07/19/19	
REV NO: 3	
SCALE N.T.S.	PAGE ELEVATED

OSHPD OPM-0486-13 DWG - 5

EQUIPMENT WEIGHTS, DIMENSIONS AND ANCHORAGE FORCES											
Model and Weight		Anchor	CG x	CG y	Eccentricities			Forces At Grade		Forces at Elevated	
Tag No.	Weight lbf	Spacing - in	in	in	ex - in	ey - in	ez - in	Vu Lbs.	Tu Lbs.	Vu Lbs.	Tu Lbs.
6-0-6 VS1	1733	26.39	19.7	17.3	-0.2	-4.8	35	845	2334	2253	4249
6-0-6 VS1 ELD	2258	26.39	17.8	18.0	-1.7	-4.1	42.0	1016	3311	3009	6707
6-0-6 VS2	2163	26.39	19.6	21.4	0.1	-0.7	35.3	933	2884	2487	5290
6-0-6 VS2 ELD	2687	26.39	18.1	18.2	-1.4	-0.9	41.1	1097	3779	3251	7739
6-0-9 VS1	2041	38.19	19.6	23.0	0.1	-5	35.7	995	2279	2653	4165
6-0-9 VS1 ELD	2566	38.19	18.0	22.4	-1.5	-5.6	41.7	1173	3042	3476	6188
6-0-9 VS2	2471	38.19	19.6	26.6	0.1	-1.4	35.8	1087	2721	2899	5014
6-0-9 VS2 ELD	2996	38.19	18.2	25.5	-1.3	-2.5	41.0	1258	3419	3728	7038
6-0-12 VS1	2372	52.16	19.6	28.4	0.1	-6.6	36.2	1186	2407	3161	4392
6-0-12 VS1 ELD	2994	52.16	17.9	28.1	-1.6	-6.8	42.4	1378	3219	4082	6545
6-0-12 VS2	2868	52.16	19.5	32.0	0	-2.9	36.3	1308	2856	3488	5267
6-0-12 VS2 ELD	3490	56.16	18.1	31.2	-1.4	-3.8	41.6	1489	3599	4411	7414

1. FORCES ARE FACTORED LOADS PER BOLT USING STRENGTH DESIGN AND INCLUDE THE FOLLOWING FACTORS: $DL = 0.9$, $S_{DS} = 2.50$, F_{pH} AT GRADE = 1.13, F_{pH} ELEVATED = 3.0 AND $F_{pV} = 0.5$ EXCEPT "ELD" MODELS AT GRADE WHERE $S_{DS} = 2.25$, $F_{pH} = 1.01$ AND $F_{pV} = 0.45$ AND (4) ANCHORS ARE USED.
2. FORCES FOR CONCRETE ANCHORAGE ARE AT STRENGTH DESIGN LEVEL AND INCLUDE A CONCRETE OVERSTRENGTH FACTOR $\Omega_0 = 1.5$ FOR SHEAR AND TENSION AT GRADE AND SHEAR FOR ELEVATED SLABS. FORCES FOR STEEL ELEMENTS ON ELEVATED FLOORS INCLUDE AN OVERSTRENGTH FACTOR $\Omega_0 = 1.0$.


BY: Sonia Eliseo

DATE: 08/09/2019

OPM-0486-13 BELIMED MST-V 6-0-X STERILIZER
DIMENSIONS AND ANCHORAGE FORCES




A Division of Tomarco Contractor Specialties
International Seismic Application Technology
 1020 Crews Road, Suite Q, Matthews, NC 28105
 704-841-4080 www.isatsb.com



DRAWN BY: WVJ
DATE: 04/23/18
REVISED BY: WVJ
DATE: 07/19/19
REV NO: 3
SCALE: N.T.S.
PAGE: DIMENSIONS


OSHPD OPM-0486-13 DWG - 6



OPM-0486-13 BELIMED MST-V 6-0-6 VS1 ELD AND
MST-V 6-0-6 VS2 ELD ANCHOR LAYOUT AT GRADE $S_{DS} = 2.50$




**A Division of Tomarco Contractor Specialties
International Seismic Application Technology**
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com



DRAWN BY: WVJ	
DATE: 04/23/18	
REVISED BY: WVJ	
DATE: 07/19/19	
REV NO: 3	
SCALE N.T.S.	PAGE V 6-0-6 ELD

OSHPD OPM-0486-13 DWG - 7




- ① BELIMED DWG NO. 000-323-947 SEISMIC BRACKET BY BELIMED (304 SS)
- ② BELIMED DWG NO. 0000-323-942 SEISMIC BRACKET BY BELIMED (304 SS)
- ③ BELIMED DWG NO. 0000-323-943 SEISMIC BRACKET BY BELIMED (304 SS)

OPM-0486-13 BELIMED MST-V 6-0-9 VS1 ELD AND
MST-V 6-0-9 VS2 ELD ANCHOR LAYOUT AT GRADE $S_{DS} = 2.50$






**A Division of Tomarco Contractor Specialties
International Seismic Application Technology**
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com




DRAWN BY: WVJ	
DATE: 04/23/18	
REVISED BY: WVJ	
DATE: 07/19/19	
REV NO: 3	
SCALE N.T.S.	PAGE V 6-0-9 ELD

OSHPD OPM-0486-13 DWG - 8



 A Division of Tomarco Contractor Specialties International Seismic Application Technology 1020 Crews Road, Suite Q, Matthews, NC 28105 704-841-4080 www.isatsb.com		 <i>WVJ 07/19/19</i>	DRAWN BY: WVJ DATE: 04/23/18
			REVISED BY: WVJ DATE: 07/19/19 REV NO: 3
SCALE N.T.S.		PAGE V 6-0-12 ELD	


OSHPD OPM-0486-13 DWG - 9



OPM-0486-13 BELIMED MST-V.6-0-X STERILIZER BRACKETS
FOR ELD MODELS AT GRADE WITH $S_{DS} = 2.50$



**A Division of Tomarco Contractor Specialties
International Seismic Application Technology**
1020 Crews Road, Suite Q, Matthews, NC 28105
704-841-4080 www.isatsb.com



DRAWN BY: WVJ
DATE: 04/23/18

REVISED BY: WVJ
DATE: 07/19/19
REV NO: 3


SCALE
N.T.S.

PAGE
BRACKETS



OSHPD OPM-0486-13 DWG - 10

EQUIPMENT WEIGHTS, DIMENSIONS AND ANCHORAGE FORCES								
ELD MODELS AT GRADE WITH $S_{DS} = 2.50$								
Model and Weight		CG x	CG y	Eccentricities			Forces At Grade	
Tag No.	Weight lbf	in	in	ex - in	ey - in	ez - in	Vu Lbs.	Tu Lbs.
6-0-6 VS1 ELD	2258	17.8	18	-1.7	-4.1	42.0	548	1575
6-0-6 VS2 ELD	2687	18	22.4	-1.4	-0.9	41.1	602	1771
6-0-9 VS1 ELD	2566	17.9	28.1	-1.5	-5.6	41.7	637	1506
6-0-9 VS2 ELD	2996	18.1	21.1	-1.3	-2.5	41.0	690	1667
6-0-12 VS1 ELD	2994	18.2	25.5	-1.6	-6.8	42.4	752	1597
6-0-12 VS2 ELD	3490	18.1	31.2	-1.4	-3.8	41.6	818	1756

- FORCES ARE FACTORED LOADS PER BOLT USING STRENGTH DESIGN AND INCLUDE THE FOLLOWING FACTORS:
DL = 0.9, $S_{DS} = 2.50$, F_{pH} AT GRADE = 1.13, F_{pH} ELEVATED = 3.0 AND $F_{pV} = 0.5$ AND (8) ANCHORS ARE USED.
- FORCES FOR CONCRETE ANCHORAGE ARE AT STRENGTH DESIGN LEVEL AND INCLUDE A CONCRETE OVERSTRENGTH FACTOR $\Omega_0 = 1.5$ FOR SHEAR AND TENSION.



OPM-0486-13 BELIMED MST-V 6-0-X ELD MODEL STERILIZER
AT GRADE DIMENSIONS AND ANCHORAGE FORCES WITH $S_{DS} = 2.50$ (8 ANCHORS)

 A Division of Tomarco Contractor Specialties International Seismic Application Technology 1020 Crews Road, Suite Q, Matthews, NC 28105 704-841-4080 www.isatsb.com		
 <i>WVJ 19 Jun 19</i>	DRAWN BY: WVJ DATE: 04/23/18	
	REVISED BY: WVJ DATE: 07/19/19 REV NO: 3	
	<table border="1"> <tr> <td>SCALE N.T.S.</td> <td>PAGE ELD GRADE</td> </tr> </table>	SCALE N.T.S.
SCALE N.T.S.	PAGE ELD GRADE	

OSHPD OPM-0486-13 DWG - 11