

Lubricant WD



Instrument lubricant formulated for automated washer/disinfectors-concentrate

Designed for automated washing of surgical instruments, other reusable items, and devices

Belimed Protect® Lubricant WD is designed to deliver superior performance in washer disinfectors. Our lubricant formula keeps instruments, hinges and box locks moving freely and provides a protective coating to help extend instrument life and prevent rusting and staining. The formula has been optimized for easy dosing and calibration in all machine types.

Belimed Protect products also have a unique color and number coded labeling system designed to reinforce what product is used at each step and reduce the risk of dosing errors. Additionally, all products have been independently tested for dermal and ocular irritation, and confirmed to be non-irritating to skin and eyes.

Product Features

- pH neutral
- Gentle on instruments
- Fragrance free
- Biodegradable and aquatic life safe

Key Formula Ingredients

- Food-grade mineral oil
- Unique surfactants for improved cleaning efficacy

Key Chemical Attributes	
pH	5.5-7.5
Density	0.995-1.005
Color	Milky

Material Compatibility

The Belimed Protect Lubricant WD demonstrates wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- Aluminum
- Anodized Aluminum
- Titanium
- Copper
- Brass
- All commonly used plastics and elastomers, including Polyurethane, Silicone and Teflon
- Latexes and Synthetics (Vinyl and Nitrile)

Storage

Do not store at extreme hot or cold temperatures

Safety Information

Before using product please refer to Instructions for Use, Safety Data Sheet and product label. Always refer to instrument IFU regarding lubricant use.

Instructions for Use

1. Automatically feed 1/12 to 1/3 oz of lubricant per gallon (0.7 to 2 ml/L) of water to the washer disinfectant
2. Follow facility policies and procedures

Product Ordering Information

	Color Code	Item Order #	Product ID (On the Label)	Pkg	Units per case
Lubricant WD	5	106-8056	2150001	1 gal	4
		106-8058	2150002	2.5 gal	2
		106-8059	2150005	5 gal	1
		106-8060	2150015	15 gal	1
		106-8051	2150030	30 gal	1
		106-8062	2150055	55 gal	1