

Enzymatic Detergent Manual Cleaning



Specially formulated multi-enzyme detergent designed to optimize the manual cleaning process – Concentrate

For manual cleaning of surgical grade instruments/reusable items, endoscopes and other devices.

Product Features:

- Proprietary, subtilisin-free formulation
- pH neutral
- Non-corrosive
- Low foaming
- Gentle on instruments
- Safe for end users
- Biodegradable and aquatic life safe

Key Formula Ingredients:

- Three (3) enzymes (protease for proteins, lipase for fatty acids and amylase for starches)
- Unique blend of surfactants solubilize bioburden
- Sequestering agent to combat hard water minerals

Key Chemical Attributes	
pH Neat	7.0-8.5
pH Use Dilution	6-8
Density	1.020-1.030
Color	Green

Material Compatibility:

The Belimed Protect® Enzymatic Detergent Manual Cleaning demonstrates a wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- Aluminum
- Anodized Aluminum
- Titanium
- Copper
- Brass
- All commonly used plastics and elastomers, including Polyurethane, Silicone, and Teflon
- Latexes and Synthetics (Vinyl and Nitrile)

Handling and Storage:

Do not store at extreme hot or cold temperatures.

Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

Instructions for Use:

1. Add 1/2 to 1 oz per gallon (4 to 8 ml/L) of water
2. For best results, ensure water temperature is 80-110°F (26-44°C). Do not exceed 140°F (60°C)
3. Soak and clean instruments according to hospital protocol
4. Rinse thoroughly

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Enzymatic Detergent Manual Cleaning Concentrate, Low Foam		106-8026	2120001	1 gal	4
		106-8028	2120002	2.5 gal	2
		106-8029	2120005	5 gal	1
		106-8030	2120015	15 gal	1
Enzymatic Detergent Manual Cleaning Concentrate, Low Foam with Fragrance		107-4535	2121001	1 gal	4
		107-4538	2121002	2.5 gal	2
		107-4539	2121005	5 gal	1

REV12.22.2020_Q0

Engineers of Confidence.