

# Enzymatic Detergent WD, Concentrate PLUS



**Specially formulated, multi-enzyme detergent designed to optimize the cleaning process in automated washer disinfectors – Concentrate**

For automated cleaning of surgical instruments, other reusable items, and devices.

**Product Features:**

- More concentrated product results in fewer changeovers and reduced storage requirements
- pH neutral
- Non-corrosive
- Gentle on instruments
- Fragrance free
- Biodegradable and aquatic life safe

**Key Formula Ingredients:**

- Three (3) Enzymes (Protease for proteins, lipase for fatty acids and amylase for starches)
- Unique blend of surfactants solubilize bioburden
- Corrosion inhibitor

Key Chemical Attributes	
pH Neat	8.0-9.0
pH 1% Dilution	6.0-8.0
Density	1.060-1.070
Color	Blue

**Material Compatibility:**

The Belimed Protect® Enzymatic Detergent WD demonstrates a wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- Aluminum
- Anodized Aluminum
- Titanium
- Copper
- Brass
- All commonly used plastics and elastomers, including Polyurethane, Silicone, and Teflon
- Latexes and Synthetics (Vinyl and Nitrile)

**Handling and Storage:**

Do not store at extreme hot or cold temperatures.

**Safety Information:**

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

**Instructions for Use:**

1. Automatically feed 1/16 to 3/8 oz of detergent per gallon (0.5 to 3 ml/L) of water
2. For best results, water temperature should be 85°F - 150°F (30°C - 65°C)
3. Follow facility policies and procedures

**Product Ordering Information:**

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Enzymatic Detergent Washer Disinfectant, Concentrate PLUS	+	106-8732	2141001	1 gal	4
		106-8734	2141002	2.5 gal	2
		106-8735	2141005	5 gal	1
		106-8736	2141015	15 gal	1



For more information, visit [belimed.com](http://belimed.com)

REV12.22.2020\_Q0

Engineers of Confidence.