


CA9618-HS1RM

RIGHT CHAMBER HORIZONTAL STERILIZER

3D IMAGE	Dwg. No. : CUT-CA9618HS1RM	Rev. : A
Typical installation of a Belimed CA9618-HS1RM Horizontal Sterilizer.	Drawn By : A. HERRON	Date : 6/23/16
	Scale : N.T.S.	Page : 1 of 11



CLEAN SIDE

- CA — COMPRESSED AIR
- CD — CONDENSATE
- CR — CHILLED WATER RETURN (OPTIONAL)
- CS — CHILLED WATER SUPPLY (OPTIONAL)
- CW — COLD WATER
- DR — FLOOR DRAIN
- DP — DATA PORT
- EL — ELECTRICITY FOR MACHINE
- SRV — SAFETY PRESSURE RELIEF VALVE
- SS — STEAM SUPPLY
- SSA — STEAM SUPPLY FROM MANIFOLD

Belimed
Infection Control
2325 Charleston Regional Parkway, Charleston, SC 29492
Tel: (800) 457-4117 Web: www.belimed.com

CA9618-HS1RM RIGHT CHAMBER HORIZONTAL STERILIZER

PLAN VIEW

Typical site preparation of a Belimed
CA9618-HS1RM Horizontal Sterilizer.


Dwg. No. :	CUT-CA9618HS1RM	Rev. :	A
Drawn By :	A. HERRON	Date :	6/23/16
Scale :	1/2" = 1'	Page :	2 of 11

A
9

PHASE 1 CONNECTION POINT FOR UTILITIES ARE TO BE PROVIDED BY OTHERS **BEFORE** EQUIPMENT INSTALLATION BEGINS. PROVIDE SHUT-OFF VALVES AT PHASE 1 CONNECTION POINTS FOR SS, CA, CW, CS, CR.

PHASE ONE
↑
PHASE TWO
↓

PHASE 2 CONNECTION OF UTILITIES TO EQUIPMENT ARE TO BE PROVIDED BY OTHERS **AFTER** EQUIPMENT HAS BEEN SET IN FINAL POSITION BY BELIMED. ALL FINAL CONNECTIONS (FROM PHASE 1 CONNECTION POINTS TO FINAL CONNECTION POINTS, PHASE 2) ARE BY OTHERS.



CA — COMPRESSED AIR
CD — CONDENSATE
CR — CHILLED WATER RETURN (OPTIONAL)
CS — CHILLED WATER SUPPLY (OPTIONAL)
CW — COLD WATER
DP — DATA PORT
DR — FLOOR DRAIN
EL — ELECTRICITY FOR MACHINE
SS — STEAM SUPPLY
SSA — STEAM SUPPLY FROM MANIFOLD

*SEE DETAIL B/9 FOR CHILLED WATER BYPASS.

LEGEND

SYMBOLS

— SERVICE SHUT-OFF VALVE, BY OTHERS
— STEAM TRAP, BY OTHERS

Belimed

Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
Tel: (800) 457-4117 Web: www.belimed.com


CA9618-HS1RM

RIGHT CHAMBER HORIZONTAL STERILIZER

FRONT ELEVATION

Typical site preparation of a Belimed
CA9618-HS1RM Horizontal Sterilizer.

Dwg. No. :	CUT-CA9618HS1RM	Rev. :	A
Drawn By :	A. HERRON	Date :	6/23/16
Scale :	1/2" = 1'	Page :	3 of 11



CA — COMPRESSED AIR
 CD — CONDENSATE
 CR — CHILLED WATER RETURN (OPTIONAL)
 CS — CHILLED WATER SUPPLY (OPTIONAL)
 CW — COLD WATER
 DP — DATA PORT
 DR — FLOOR DRAIN
 EL — ELECTRICITY FOR MACHINE
 SS — STEAM SUPPLY
 SSA — STEAM SUPPLY FROM MANIFOLD

LEGEND

SYMBOLS

— SERVICE SHUT-OFF VALVE, BY OTHERS
 — STEAM TRAP, BY OTHERS

Belimed

Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
 Tel: (800) 457-4117 Web: www.belimed.com

CA9618-HS1RM

RIGHT CHAMBER HORIZONTAL STERILIZER

SIDE ELEVATION

Typical site preparation of a Belimed
 CA9618-HS1RM Horizontal Sterilizer.

Dwg. No. : CUT-CA9618HS1RM


Rev. : A

Drawn By : A. HERRON

Date : 6/23/16

Scale : 1/2" = 1'

Page : 4 of 11



CA — COMPRESSED AIR
 CD — CONDENSATE
 CR — CHILLED WATER RETURN (OPTIONAL)
 CS — CHILLED WATER SUPPLY (OPTIONAL)
 CW — COLD WATER
 EL — ELECTRICITY FOR MACHINE
 SRV — SAFETY PRESSURE RELIEF VALVE
 SS — STEAM SUPPLY
 SSA — STEAM SUPPLY FROM MANIFOLD

LEGEND

SYMBOLS

— SERVICE SHUT-OFF VALVE, BY OTHERS
 — STEAM TRAP, BY OTHERS

Belimed

Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
 Tel: (800) 457-4117 Web: www.belimed.com


CA9618-HS1RM

RIGHT CHAMBER HORIZONTAL STERILIZER


TYPICAL ISOMETRIC VIEW

Typical utility piping diagram of a Belimed CA9618-HS1RM Horizontal Sterilizer.

Dwg. No. :	CUT-CA9618HS1RM	Rev. :	A
Drawn By :	A. HERRON	Date :	6/23/16
Scale :	N.T.S.	Page :	5 of 11



CLEAN SIDE



STERILIZER TRANSPORT CART AND LOADING
RACKS SOLD SEPARATELY. SEE CART CUT
SHEET FOR ADDITIONAL DIMENSIONS.

Belimed

Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
Tel: (800) 457-4117 Web: www.belimed.com


CA9618-HS1RM

RIGHT CHAMBER HORIZONTAL STERILIZER

PLAN VIEW


Typical equipment layout of a Belimed
CA9618-HS1RM Horizontal Sterilizer.

Dwg. No. :	CUT-CA9618HS1RM	Rev. :	A
Drawn By :	A. HERRON	Date :	6/23/16
Scale :	1/2" = 1'	Page :	6 of 11



LEGEND

HATCH TYPES

 STAINLESS STEEL TRIM, BY BELIMED

Belimed

Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
Tel: (800) 457-4117 Web: www.belimed.com

CA9618-HS1RM

RIGHT CHAMBER HORIZONTAL STERILIZER


FRONT ELEVATION

Typical equipment layout of a Belimed CA9618-HS1RM Horizontal Sterilizer.

Dwg. No.: CUT-CA9618HS1RM Rev.: A

Drawn By: A. HERRON Date: 6/23/16


Scale: 1/2" = 1' Page: 7 of 11



CA9618-HS1RM


RIGHT CHAMBER HORIZONTAL STERILIZER

SIDE ELEVATION	Dwg. No. : CUT-CA9618HS1RM	Rev. : A
Typical equipment layout of a Belimed CA9618-HS1RM Horizontal Sterilizer.	Drawn By : A. HERRON	Date : 6/23/16
	Scale : 1/2" = 1'	Page : 8 of 11



NOTE: REFER TO LOCAL BUILDING CODES AND/OR SEISMIC ANCHORAGE REQUIREMENTS FOR FINAL ATTACHMENT METHOD.

1. UTILITY ATTACHMENT ASSEMBLY SHALL FALL WITH IN THE FOOTPRINT OF THE SPECIFIED EQUIPMENT.
2. REFER TO EQUIPMENT PLUMBING, MECHANICAL AND ELECTRICAL PLANS FOR UTILITY LOCATIONS.
3. SEE TECHNICAL SPECIFICATIONS FOR UTILITY SIZING, INPUT REQUIREMENTS, AND FINAL CONNECTIONS.




MANUAL SERVICE VALVES, BY OTHERS.


CW - COLD WATER
CS - CHILLED WATER SUPPLY
CR - CHILLED WATER RETURN

A RECOMMENDED TYPICAL CEILING HUNG UTILITY ATTACHMENT DETAIL (SIDE VIEW)
N.T.S.

B TYP. (OPTIONAL) BY PASS DETAIL CHILLED WATER CONNECTION
N.T.S.




(TO BE INSTALLED BY OTHERS WITH 1/2\"/>



C TYPICAL STEAM MANIFOLD DETAIL
N.T.S.

D TYP. TRIM ATTACHMENT DETAIL HORIZONTAL STERILIZER
N.T.S.



E BASE COPING AT WALL OPENING
N.T.S.

CA9618-HS1 STERILIZER - TECHNICAL SPECIFICATIONS

REQUIRED INPUTS

Symbol	Description	Value	Utility Notes	Connection Method	Final Utility Conn. By:	
					Others	Belimed
SS	Utility:	Steam	Steam must be dry, saturated, and free from all impurities such as gases, oil, dirt, or other contaminants. Provide adequate steam strainer typical for sterilizers. All steam drops must be piped from top of steam supply line.	Hard Pipe	X	
	Connection Size/Type:	1 1/2 inch / MNPT				
	Pressure:	Min: 40 Max: 50 psig				
	Peak Consumption:	6.25 lbs/min				
CW	*Consumption per Load:	88 lbs	Cold water is used to prime the condensate pump and cool the condensate. Colder water temperatures will result in less water consumption. Optimal temperature is 60°F. Temperature must be less than or equal to 85°F. TDS must be less than or equal to less than 400 ppm. Back flow prevention must be supplied by customer, as per local code.	Flex		X
	Utility:	Cold Water				
	Connection Size/Type:	1 inch / MNPT				
	Pressure:	Min: 30 Max: 70 psig				
CA	Peak Consumption:	14.0 gal/min	Air must be dry and oil free. Used to operate pneumatic control system.	Flex		X
	*Consumption per Load:	167 gal				
	Utility:	Compressed Air				
	Connection Size/Type:	1/4 inch / MNPT				
EL	Pressure:	Min: 70 Max: 100 psig	Junction Box with Receptacle. NEMA 15-30R to NEC local codes. Grainger part No. 5Z889.	Receptacle		X
	Peak Consumption:	6 CFM				
	*Consumption per Load:	14 ft³				
	Utility:	Electric				
SRV1	Connection Type:	Receptacle	Provide sufficient air exhaust to remove 4 kW heat dissipation from service room.		X	
	Voltage (Nominal):	208 VAC				
	Voltage (Phase):	3 Phase				
	Voltage (Frequency):	60 Hz				
SRV2	Circuit Protection:	30 Amps	Provide connection to safety relief valve and pipe to drain, atmosphere, or as required by local code.	Hard Pipe with Union	X	
	Peak Consumption:	4.9 kW				
	*Consumption per Load:	4.2 kW-h				
	Utility:	Equipment Area Exhaust				
SSA	Connection Type:	Exhaust return grille	Pipe from manifold to sterilizers by others.	Hard Pipe	X	
	Utility (Chamber):	Safety Pressure Relief Valve				
	Outlet size:	1 1/2 inch / MNPT				
	Utility (Jacket):	Safety Pressure Relief Valve				
DP	Outlet size:	1 inch / MNPT	Location: Ethernet port should be located within 5 ft of machine within the enclosed service area of the sterilizer.		X	
	Utility:	Data Port				
	Connection Method:	TCP/IP Ethernet to hospital network				
	Connection Quantity:	One (1) RJ45 Ethernet port per sterilizer				
CD	Facility Wiring:	*Cat-5, Cat-6 or Cat-6e wire	* (By others)	Hard Pipe	X	
	NIC/Com Speed:	10/100 Mbps Networking Port (RJ45)				
	Broadcast Method:	Unicast				
	MTU Size:	1500				
DR	Final Connection:	Belimed provided patch cables for final connection	Each RJ45 Ethernet Port with individual static IP address		X	
	Port IP Setup:	Each RJ45 Ethernet Port with individual static IP address				
	Utility:	Open Drain				
	Connection Size:	2 inch				
CD	Discharge Flow Rate:	16.4 gal/min	Condensate removal by customer per local code.	Hard Pipe	X	
	*Consumption per Load:	177.6 gal (w/o chilled water)				
	Medium:	Noise Radiation				
	Radiation (Net):	57 dBA				
CD	Medium:	Heat Radiation	Heat dissipated into service room with ambient temperature of 86°F.			
	Radiation (Net):	4 kWh				
	Utility:	Condensate				
	Connection Size/Type:	3/4 inch / FNPT				

ESTIMATED EMISSIONS

DR	Utility:	Open Drain	Provide flush mounted floor sink and trap per local code.	2" Copper Pipe	X	
	Connection Size:	2 inch				
	Discharge Flow Rate:	16.4 gal/min				
	*Consumption per Load:	177.6 gal (w/o chilled water)				
CD	Medium:	Noise Radiation	In accordance with ISO3746: 1995 Testing Standards.			
	Radiation (Net):	57 dBA				
	Medium:	Heat Radiation				
	Radiation (Net):	4 kWh				
CD	Utility:	Condensate	Condensate removal by customer per local code.	Hard Pipe	X	
	Connection Size/Type:	3/4 inch / FNPT				

CONTINUED ON PAGE 11

Belimed

Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
Tel: (800) 457-4117 Web: www.belimed.com

CA9618-HS1RM RIGHT CHAMBER HORIZONTAL STERILIZER

TECHNICAL SPECIFICATIONS

Technical Data for the installation of a Belimed CA9618-HS1RM Horizontal Sterilizer.

Dwg. No. :	CUT-CA9618HS1RM	Rev. :	A
Drawn By :	A. HERRON	Date :	6/23/16
Scale :	NONE	Page :	10 of 11

Optional utilities for chilled water cooling (used for reduction of cold water consumption).

CS	Utility:	Chilled Water Supply	Feed water at ~50°F. Approximate $\Delta T = 15/25$ K. Cooling Capacity: average/peak	Insulated Hard Pipe	X	
	Connection Size/Type:	1 1/4 inch / MNPT				
	Pressure:	Min: 30 Max: 70 psig				
	Peak Flow Rate:	19.8 gal/min				
	Peak Consumption:	678 kBTU/hr				
	*Consumption per Load:	169 kBTU				
CR	Utility:	Chilled Water Return	Pressure difference ΔP CS/CR must be greater than 7 PSIG.	Insulated Hard Pipe	X	
	Connection Size/Type:	1 1/4 inch / MNPT				
CW	Utility:	Cold Water	Cold water is used to prime the condensate pump. Water should be softened to 0-60 ppm. Back flow prevention must be supplied by customer, as per local code.	Flex		X
	Connection Size/Type:	1/2 inch / MNPT				
	Peak Consumption:	2.2 gal/min				
	*Consumption per Load:	2.6 gal				

PHYSICAL SPECIFICATIONS

	Parameter:	Carry-In Dimensions			Machine Alone (i.e. removed from packaging)
		Chamber	Tech Frame	Face Panel	
	Height:	70.4 in	72 in	78.7 in	
	Width:	39 in	31 in	70.9 in	
	Length:	90 in	40 in	31.5 in	
	Parameter:	Chamber Usable Dimensions			
	Height:	42.5 in			
	Width:	26 in			
	Length:	79.25 in			
	Parameter:	Weights			
	Transport Weight:	5035 lbs			
	Operating Weight:	5236 lbs			
	OSHPD Number:	Pending Approval			

*Consumptions calculated under pre-vac cycle with 43 psig at 270° F and water temperature of 60° F.

Additional Information:

- MNPT = Male pipe thread - FNPT = Female pipe thread
- Others must provide service valves and strainers for steam and water supply.
- All utility supplies are other's responsibility.
- Left or Right side chamber is available.

TS-CA9618-HS1_RevA



Infection Control

2325 Charleston Regional Parkway, Charleston, SC 29492
Tel: (800) 457-4117 Web: www.belimed.com

CA9618-HS1RM
RIGHT CHAMBER HORIZONTAL STERILIZER

TECHNICAL SPECIFICATIONS CONT Technical Data for the installation of a Belimed CA9618-HS1RM Horizontal Sterilizer.	Dwg. No. :	CUT-CA9618HS1RM	Rev. :	A
	Drawn By :	A. HERRON	Date :	6/23/16
	Scale :	NONE	Page :	11 of 11