



OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

APPLICATION FOR OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM)

OFFICE USE ONLY
APPLICATION #: OPM-0219-13

OSHPD Preapproval of Manufacturer's Certification (OPM)

Type: [X] New [ ] Renewal [ ] Update to Pre-CBC 2013 OPA Number:

Manufacturer Information

Manufacturer: Belimed, Inc.

Manufacturer's Technical Representative: Nelson Garrido

Mailing Address:

Telephone: 843-216-7424 Email: nelson.garrido@belimed.com

Product Information

Product Name: Ultrasonic Washer

Product Type: Laboratory Equipment OPM-0219-13

Product Model Number: 200c, 400, & UW565

General Description: Floor mounted ultrasonic washers with seismic mounting brackets.

Applicant Information

Applicant Company Name: Belimed, Inc.

Contact Person: Nelson Garrido

Mailing Address: 2325 Charleston Regional Pkwy., Charleston, SC. 29492

Telephone: 843-216-7424 Email: nelson.garrido@belimed.com

I hereby agree to reimburse the Office of Statewide Health Planning and Development review fees in accordance with the California Administrative Code, 2013.

Signature of Applicant: [Signature] Date: 4/1/2015
Manager of Engineering, Quality and Regulatory Affairs
Company Name: Belimed, Inc.

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"





OFFICE OF STATEWIDE HEALTH PLANNING AND DEVELOPMENT
FACILITIES DEVELOPMENT DIVISION

Registered Design Professional Preparing Engineering Recommendations

Company Name: Makeitright Inc.

Name: Joseph La Brie California License Number: S3566

Mailing Address: 55 E. Huntington Drive Ste 277, Arcadia, CA 91006

Telephone: (626) 445-0366 Email: labrie@makeitright.net

OSHPD Special Seismic Certification Preapproval (OSP)

- Special Seismic Certification is preapproved under OSP- (Separate application for OSP is required)
Special Seismic Certification is not preapproved

Certification Method(s)

- Testing in accordance with: ICC-ES AC156 FM 1950-10
Other\* (Please Specify): Ultrasonic washer does not require special seismic certification.

\*Use of criteria other than those adopted by the California Building Standards Code, 2013 (CBSC 2013) for component supports and attachments are not permitted. For distribution system, interior partition wall, and suspended ceiling seismic bracings, test criteria other than those adopted in the CBSC 2013 may be used when approved by OSHPD prior to testing.

- Analysis
Experience Data
Combination of Testing, Analysis, and/or Experience Data (Please Specify):

List of Attachments Supporting the Manufacturer's Certification

- Test Report Drawings Calculations Manufacturer's Catalog
Other(s) (Please Specify):


OFFICE USE ONLY - OSHPD APPROVAL VALID FOR CBC 2013 ONLY

Signature: [Signature] Date: 06/25/2015
Print Name: William Staehlin
Title: SSE
Condition of Approval (if applicable):



# Belimed

Infection Control



## ULTRASONIC WASHER



MODELS: 200c, 400, & UW565

OPM-0219-13

TABLE OF CONTENTS

BELIMED, INC. – ULTRASONIC WASHER  
MODEL NOS. 200c, 400, & UW565



TC-1	TABLE OF CONTENTS
GN-1	GENERAL NOTES
GN-2	GENERAL NOTES
1.0	MODEL 200c DETAIL
1.1	MODEL 400 DETAIL
1.2	MODEL UW565 DETAIL
2.0	SEISMIC MOUNT DETAIL
3.0	ANCHORAGE TO CONCRETE SLAB ON GRADE DETAIL
4.0	ANCHORAGE TO TYP. "W" DECK DETAIL
4.1	ANCHORAGE TO TYP. "B" DECK DETAIL



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:  <b>JOSEPH L. LA BRIE</b> Structural Engineer No. SE 3566   hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947	MANUFACTURER:   Infection Control	PRODUCT: ULTRASONIC WASHER Models.: 200c, 400, & UW565	DATE: 04/01/15
		PAGE TITLE: TABLE OF CONTENTS	PAGE: TC-1

GENERAL NOTES

1. THIS OSHPD PREAPPROVAL OF MANUFACTURER'S CERTIFICATION (OPM) IS BASED ON THE 2013 CALIFORNIA BUILDING CODE (CBC). THE DEMAND (DESIGN FORCES) FOR USE WITH THIS OPM SHALL BE BASED ON THE 2013 CALIFORNIA BUILDING CODE (CBC).
2. SITE VERIFICATION IS REQUIRED. THE CONTRACTOR SHALL VERIFY ALL EXISTING CONDITIONS IN THE AREA(S) OF WORK PRIOR TO THE BEGINNING OF THE PROJECT. DO NOT SCALE THE DRAWINGS; ALL DIMENSIONS MUST BE VERIFIED IN THE FIELD. THE STRUCTURAL ENGINEER SHALL BE NOTIFIED, IMMEDIATELY, IF ANY DISCREPANCIES ARE FOUND.
3. DESIGN CRITERIA:
  - a.) SLAB ON GRADE  
 $S_{DS} = 2.5 \text{ MAX.}; l = 1.5 ; z/h = 0$  (GROUND OR SLAB ON GRADE)  
UPPER FLOOR  
 $S_{DS} = 2.5 \text{ MAX.}; l = 1.5 ; z/h \leq 1$  (UPPER FLOOR)
  - b.) PER ASCE 7-10 INCLUDING SUPPLEMENT NO. 1 & TABLE 13.6-1 :  
 $\alpha_p = 1.0 ; R_p = 1.5 ; \Omega_o = 1.5$  (APPLY  $\Omega_o$  FACTOR FOR ANCHORAGE TO CONCRETE)
4. CENTER OF GRAVITY (C.G.) WEIGHT IS A MAXIMUM.
5. STRUCTURAL ENGINEER OF RECORD FOR THE BUILDING SHALL :
  - a. CHECK THAT THE EXISTING/NEW STRUCTURE WILL BE ABLE TO SUPPORT THE MAXIMUM WEIGHTS/FORCES SHOWN IN ADDITION TO ANY OTHER LOADS TO THE STRUCTURE. PROVIDE STRENGTHENING OF STRUCTURE AS REQUIRED.
  - b. CHECK THAT THE FLOOR OR DECK ANCHORS ARE LOCATED AT AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS. ALL MINIMUM EDGE DISTANCE AND SPACING SHALL BE IN ACCORDANCE WITH THE REQUIREMENTS GIVEN IN THE ICC REPORT FOR THE SPECIFIC ANCHORS USED ON THIS OPM. (SEE TABLE ON NEXT PAGE FOR ANCHOR MINIMUM SPACING & EDGE DISTANCE REQUIREMENTS)
  - c. CHECK THAT THE INSTALLATION, SUPPORT AND ATTACHMENTS OF THE UNIT COMPLIES WITH THE 2013 CALIFORNIA BUILDING CODE AND WITH THE DETAILS SHOWN IN THIS PREAPPROVAL.
  - d. VERIFY THAT THE ACTUAL EQUIPMENT'S MODEL NUMBER, OPERATING WEIGHT, CENTER OF GRAVITY (C.G.) LOCATION, ANCHOR LAYOUT, MATERIAL & ASTM GRADE OF THE EQUIPMENT IS THE SAME SHOWN ON THIS OPM PREAPPROVAL.
  - e. VERIFY THAT THE PROJECT SPECIFIC  $S_{DS}$  AND  $z/h$  VALUES RESULT IN SEISMIC FORCES DOES NOT EXCEED THE VALUES SHOWN ON THIS OPM.



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:  <b>JOSEPH L. LA BRIE</b> Structural Engineer No. SE 3566   hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947	MANUFACTURER:   Infection Control	PRODUCT: ULTRASONIC WASHER Models.: 200c, 400, & UW565	DATE: 04/01/15
		PAGE TITLE: GENERAL NOTES	PAGE: GN-1

GENERAL NOTES (CONTINUED)

6. EXPANSION ANCHORS SHALL BE HILTI KB-TZ ANCHORS (PER ICC-ESR-1917 DATED AUGUST 2014) AND INSTALLED IN NORMAL WEIGHT CONCRETE & SAND LIGHT WEIGHT CONCRETE (AT UNDERSIDE OF DECK) WHERE OCCURS. CARBON STEEL FOR INDOOR APPLICATIONS. MINIMUM EMBEDMENT OF ALL BOLTS AND TEST LOADS (UNLESS NOTED OTHERWISE ON DETAIL):

TEST LOADS FOR HILTI KB-TZ EXPANSION ANCHORS  
(f'c = 3000 PSI MIN.) INSTALLED TO N.W. CONCRETE SLAB & SAND L.W. CONCRETE AS OCCURS

BOLT SIZE	APPLICATION	MIN. EFFECTIVE EMBEDMENT	MIN. ANCHOR SPACING	MIN. EDGE DISTANCE	INSTALLATION TORQUE (LBS.-FT.)
3/8"φ	UNDERSIDE OF DECK	2"	6"	6"	25
1/2"φ	TOP OF SLAB ON GRADE	2"	6"	6"	40

- a. WHEN INSTALLING DRILLED-IN ANCHORS IN EXISTING NON-PRESTRESSED REINFORCED CONCRETE, USE CARE AND CAUTION TO AVOID CUTTING OR DAMAGING THE THE EXISTING REINFORCING BARS. MAINTAIN A MINIMUM CLEARANCE OF ONE INCH BETWEEN THE REINFORCEMENT AND THE DRILLED-IN ANCHOR.
- b. AFTER A MIN. OF (24) HRS. HAVE ELAPSED, ALL POST INSTALLED ANCHORS (LOADED IN EITHER PULL OUT OR SHEAR) SHALL BE TORQUED OR TENSION TESTED. WHEN POST-INSTALLED ANCHORS ARE USED FOR NON STRUCTURAL APPLICATIONS SUCH AS EQUIPMENT ANCHORAGE, 50% OR ALTERNATE BOLTS IN A GROUP, INCLUDING AT LEAST ONE-HALF THE ANCHORS IN EACH GROUP SHALL BE TESTED. IF THERE ARE ANY FAILURES, THE IMMEDIATELY ADJACENT BOLTS MUST THEN ALSO BE TESTED.
- c. THE FOLLOWING CRITERIA APPLY FOR THE ACCEPTANCE OF INSTALLED ANCHORS:


HYDRAULIC RAM METHOD: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE APPLICABLE TEST LOAD. FOR WEDGE AND SLEEVE TYPE ANCHORS, A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER UNDER THE NUT BECOMES LOOSE. FOR ADHESIVE ANCHORS, WHERE OTHER THAN BOND IS BEING TESTED, THE DEVICE SHALL NOT RESTRICT THE CONCRETE SHEAR CONE TYPE FAILURE MECHANISM FROM OCCURRING.

TORQUE WRENCH METHOD: THE APPLICABLE TEST TORQUE MUST BE REACHED WITHIN THE FOLLOWING LIMITS:  
WEDGE OR SLEEVE TYPE: ONE-HALF (1/2) TURN OF THE NUT

- d. OWNER'S REPRESENTATIVE IS RESPONSIBLE FOR ALL ANCHOR TESTING.
- e. ALL TESTS SHALL BE PERFORMED IN THE PRESENCE OF THE INSPECTOR OF RECORD AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO THE ENFORCEMENT AGENCY.

7. BOLTS THROUGH CONCRETE ON METAL DECK:

- a. BOLTS SHALL BE TORQUED BY 3/4 TURN OF THE NUTS AFTER THE SNUG TIGHT (THE NUG-TIGHT CONDITION IS DEFINED AS THE TIGHTNESS REQUIRED TO BRING THE CONNECTED PLIES INTO FIRM CONTACT) CONDITION IS ACHIEVED, UNLESS OTHERWISE NOTED.
- b. THROUGH BOLT HOLES SHALL BE 1/16" LARGER THAN BOLT SIZE (HOLE SIZE = BOLT SIZE + 1/16) FOR CONCRETE
- c. THROUGH-BOLTS IN CONCRETE SHALL RECEIVE SPECIAL INSPECTION AND TESTING (THROUGH BOLTS WITH STEEL TO STEEL CONNECTION IN TENSION DO NOT REQUIRE TENSION TESTING) IN ACCORDANCE WITH REQUIREMENTS FOR POST-INSTALLED ANCHORS.



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:

JOSEPH L. LA BRIE  
Structural Engineer  
No. SE 3566



MANUFACTURER:



PRODUCT:

ULTRASONIC WASHER  
Models.: 200c, 400, & UW565

PAGE TITLE:

GENERAL NOTES

DATE:

04/01/15


PAGE:

GN-2

**MODEL 200c:**


**NOTES:**

1. OPERATIONAL WEIGHT OF MODEL 200c = 470 LBS. MAX.
2. EXTERIOR SHELL OF UNIT IS 16 GA. 304 STAINLESS STEEL.



**FRONT VIEW**



**RIGHT SIDE VIEW**



**PLAN VIEW @ BASE**

CONC. SLAB OR CONC./METAL DECK;  
SEE PAGE 3.0 OR 4.0 AS APPLICABLE

SEISMIC MOUNT, TYP. Ⓞ (4) PLACES,  
(SEE 3.0, 4.0 & 4.1)



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:

**JOSEPH L. LA BRIE**  
Structural Engineer  
No. SE 3566



MANUFACTURER:



PRODUCT:

ULTRASONIC WASHER  
Models.: 200c, 400, & UW565

PAGE TITLE:

MODEL 200c DETAIL


DATE:

04/01/15

PAGE:

1.0


**MODEL 400:**



- NOTES:**
1. OPERATIONAL WEIGHT OF MODEL 400 = 745 LBS. MAX.
  2. EXTERIOR SHELL OF UNIT IS 16 GA. 304 STAINLESS STEEL.

**FRONT VIEW**


**RIGHT SIDE VIEW**



**PLAN VIEW @ BASE**

SEISMIC MOUNT, TYP. @ (4) PLACES, (SEE 3.0, 4.0 & 4.1)

CONC. SLAB OR CONC./METAL DECK; SEE PAGE 3.0 OR 4.0 AS APPLICABLE






DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:  
**JOSEPH L. LA BRIE**  
 Structural Engineer  
 No. SE 3566  
  
 hospital building design professionals  
 55 E. Huntington Dr., Suite 277  
 Arcadia, CA 91006  
 PHN: (626) 445-0366  
 FAX: (626) 445-2947

MANUFACTURER:  
  
 Infection Control


PRODUCT:  
 ULTRASONIC WASHER  
 Models.: 200c, 400, & UW565  
 PAGE TITLE:  
 MODEL 400 DETAIL

DATE:  
 04/01/15  
 PAGE:  
**1.1**


**MODEL UW565:**



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:  <b>JOSEPH L. LA BRIE</b> Structural Engineer No. SE 3566   hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947	MANUFACTURER:   Infection Control	PRODUCT: ULTRASONIC WASHER Models.: 200c, 400, & UW565	DATE: 04/01/15
		PAGE TITLE: MODEL UW565 DETAIL	PAGE: 1.2




**PLAN VIEW**



**FRONT VIEW**

**SIDE VIEW**



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:

**JOSEPH L. LA BRIE**  
Structural Engineer  
No. SE 3566



**makeitright™**  
hospital building design professionals  
55 E. Huntington Dr., Suite 277  
Arcadia, CA 91006  
PHN: (626) 445-0366  
FAX: (626) 445-2947

MANUFACTURER:



**Belimed**  
Infection Control

PRODUCT:

ULTRASONIC WASHER  
Models.: 200c, 400, & UW565

PAGE TITLE:

SEISMIC MOUNT DETAIL


DATE:

04/01/15

PAGE:


2.0

**SLAB ON GRADE ANCHORAGE (z/h = 0):**




**MAX. LOADS ON ATTACHMENT:**

$\Omega_0 T_u = 681 \text{ LBS.}$   
 $\Omega_0 V_u = 469 \text{ LBS}$



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:  <b>JOSEPH L. LA BRIE</b> Structural Engineer No. SE 3566   hospital building design professionals 55 E. Huntington Dr., Suite 277 Arcadia, CA 91006 PHN: (626) 445-0366 FAX: (626) 445-2947	MANUFACTURER:   Infection Control	PRODUCT: ULTRASONIC WASHER Models.: 200c, 400, & UW565	DATE: 04/01/15
		PAGE TITLE: ANCHORAGE TO CONCRETE SLAB ON GRADE DETAIL	PAGE: 3.0

UPPER FLOOR ANCHORAGE ON TYP. "W" DECK (z/h <= 1.0):




NOTES:

1. THIS DETAIL CAN ALSO BE USED ON MIN. 4" THICK N.W. CONCRETE FLAT SLAB FLOORS ( $f'_c = 3000$  psi MIN.).

MAX. LOADS ON ATTACHMENT:

$T_u = 1270$  LBS.  
 $\Omega_o V_u = 1250$  LBS.



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:

JOSEPH L. LA BRIE  
 Structural Engineer  
 No. SE 3566



MANUFACTURER:



PRODUCT:

ULTRASONIC WASHER  
 Models.: 200c, 400, & UW565

PAGE TITLE:

ANCHORAGE TO TYP. "W"  
 DECK DETAIL



DATE:

04/01/15

PAGE:

4.0

UPPER FLOOR ANCHORAGE ON TYP. "B" DECK (z/h <= 1.0):



DESIGN PROFESSIONAL PREPARING ENGINEERING RECOMMENDATIONS:  <b>JOSEPH L. LA BRIE</b> Structural Engineer No. SE 3566   <small>hospital building design professionals          55 E. Huntington Dr., Suite 277          Arcadia, CA 91006          PHN: (626) 445-0366          FAX: (626) 445-2947</small>	MANUFACTURER:   <b>Belimed</b> Infection Control	PRODUCT: ULTRASONIC WASHER Models.: 200c, 400, & UW565	DATE: 04/01/15
		PAGE TITLE: ANCHORAGE TO TYP. "B" DECK DETAIL	PAGE: <b>4.1</b>