



**Office of Statewide Health Planning and Development**  
**ANCHORAGE PRE-APPROVAL**

**OPA-2509-10**

Equipment Manufacturer: Belimed, Inc.

Equipment Type: 400 Series Ultrasonic Cleaner

**GENERAL NOTES**


**1. EXPANSION ANCHORS:**

- (a) ATTACHMENT IS TO BE MADE WITH THE ANCHORS LISTED BELOW AND INSTALLED AS DESCRIBED IN THE CORRESPONDING ICC REPORT.

| Anchor Diameter | Concrete Type | Min. f'c (psi) | Anchor Type        | ICC Report No. | Min. Embed. | Min. Spacing | Min. Edge Dist. | Min. Conc. Thickness | Installation. Torque | Test Loads                     |
|-----------------|---------------|----------------|--------------------|----------------|-------------|--------------|-----------------|----------------------|----------------------|--------------------------------|
| 3/8"            | Hardrock      | 2500           | Hilti Kwik Bolt TZ | ESR-1917       | 2"          | 6"           | 9"              | 4"                   | 25 Ft-Lbs            | Direct Pull Tension - 1580 lbs |

**2. TESTING OF EXPANSION ANCHORS PER 2010 CBC, 1916A-8 : TENSION TESTING SHALL BE DONE IN THE PRESENCE OF THE SPECIAL INSPECTOR AND A REPORT OF THE TEST RESULTS SHALL BE SUBMITTED TO OSHPD**

- (a) AFTER AT LEAST 24 HOURS HAVE ELAPSED SINCE INSTALLATION, DIRECT PULL TENSION TEST AT LEAST 50% OF THE ANCHORS.
- (b) ACCEPTANCE CRITERIA: THE ANCHOR SHOULD HAVE NO OBSERVABLE MOVEMENT AT THE TEST LOAD. A PRACTICAL WAY TO DETERMINE OBSERVABLE MOVEMENT IS THAT THE WASHER BECOMES LOOSE.
- (c) IF ANY ANCHOR FAILS, TEST ALL ANCHORS.





**EASE**  
EQUIPMENT ANCHORAGE & SEISMIC ENGINEERING

2801 Connery Way, Suite B

Missoula, MT 59808

Phn: (406) 541-EASE (3273) Fax: (406) 541-3274

Sheet 2 of 4

## Office of Statewide Health Planning and Development

### ANCHORAGE PRE-APPROVAL

# OPA-2509-10

Equipment Manufacturer: Belimed, Inc.

Equipment Type: 400 Series Ultrasonic Cleaner


#### GENERAL NOTES (CONTINUED)

3. FORCES PER ASCE 7-05 SECTION 13.3.1, EQUATIONS 13.3-1, 13.3-2 & 13.3-3, WHERE  $S_D = 1.93$ ,  $a_p = 1.0$ ,  $I_p = 1.5$  &  $R_p = 2.5$
4. THIS PRE-APPROVAL CONFORMS TO THE 2010 CALIFORNIA BUILDING CODE.
5. THE DETAILS IN THIS PRE-APPROVAL MAY BE USED AT ANY LOCATION IN THE STATE OF CALIFORNIA.  
THE ELEVATED FLOOR DETAILS MAY BE USED AT ANY HEIGHT IN A BUILDING.
6. ALL ANCHOR FORCES SHOWN ON THE DRAWINGS ARE WORKING LOADS (AS OPPOSED TO STRENGTH LEVEL LOADS)  
AND MAY BE USED FOR ALLOWABLE STRESS DESIGN.
7. THIS PRE-APPROVAL COVERS THE ANCHORAGE OF THE UNIT ONLY.

#### RESPONSIBILITIES OF THE STRUCTURAL ENGINEER OF RECORD

8. VERIFY THAT THE CONCRETE SLAB WHICH THE EQUIPMENT IS ANCHORED TO IS NOT CRACKED.
9. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY SLAB EDGES OR OPENINGS.
10. VERIFY THAT THE ANCHORS ARE AN ADEQUATE DISTANCE FROM ANY NEW OR EXISTING ANCHORS.
11. DESIGN BACKING BARS, STUDS, ETC. WHICH THE UNITS ARE ATTACHED TO AS NOTED ON THE DRAWINGS.  
THE SEOR SHALL ALSO VERIFY THE ADEQUACY OF THE STRUCTURES (SUCH AS WALLS AND FLOORS)  
WHICH SUPPORT THE UNITS FOR THE LOADS IMPOSED ON THEM BY THE UNITS AS WELL AS ALL OTHER LOADS.
12. PROVIDE ANY SUPPORTING STRUCTURE REQUIRED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN ADDITION TO ALL OTHER LOADS.
13. VERIFY THAT THE INSTALLATION IS IN CONFORMANCE WITH THE 2010 CBC AND WITH THE DETAILS SHOWN IN THIS PRE-APPROVAL.  
VERIFY THAT THE ACTUAL EQUIPMENT'S WEIGHT, CG LOCATION, ANCHOR LOCATIONS, ANCHOR DETAILS AND THE MATERIAL AND GAGE OF THE UNIT WHERE ATTACHMENTS ARE MADE AGREE WITH THE INFORMATION SHOWN ON THE PRE-APPROVAL DOCUMENTS.



|   |                    |
|---|--------------------|
| <b>A P P R O V E D</b>  |                    |
| Fixed Equipment Anchorage<br>Office of Statewide Health Planning and Development      |                    |
|  | <b>OPA-2509-10</b> |
| Pre-approval Program Manager:<br>Anthony R. Pike<br>(916) 440-8470                    |                    |
| <b>P E N D I N G</b>  |                    |
| Reviewed By: Anthony R. Pike  | 4/7/10             |

**BELIMED, INC.**

**400 SERIES ULTRASONIC CLEANER**

DES. **R. LA BRIE**

JOB NO. **11-1018**

DATE **4/7/10**


SHEET

**3**

OF **4** SHEETS


SEISMIC ANCHORAGE

SLAB ON GRADE



C.G. WT. = 145 LBS.  
(INCLUDING CONTENTS)

(12 GA., 50 ksi (MIN))  
BASE FRAME  
w/ 4- 3/8"φ HILTI KB-TZ  
EXPANSION ANCHORS  
(MIN. EMBED. (het) = 2")



FRONT ELEVATION

SIDE ELEVATION

NOTES:

1. ANCHORAGE DESIGN PER 2010 CALIFORNIA BUILDING CODE - SECTION 1613A AND ASCE 7-05 SECTIONS 12 AND 13. ALLOWABLE STRESS DESIGN IS USED.

HORIZONTAL FORCE ( $E_h$ ) =  $0.61 W_p$  ( $S_{Ds} = 1.93$ ,  $a_p = 1.0$ ,  $I_p = 1.5$ ,  $R_p = 2.5$ )

VERTICAL FORCE ( $E_v$ ) =  $0.27 W_p$

$T_{MAX} = 149$  LBS/BOLT

$V_{MAX} = 151$  LBS/BOLT

2. CENTER OF GRAVITY (C.G.) WEIGHT IS A MAXIMUM. THIS PRE-APPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
3. ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN ADDITION TO ALL OTHER LOADS.
4. SEE GENERAL NOTES: SHEETS 1-2



**BELIMED, INC.**

**400 SERIES ULTRASONIC CLEANER**

DES. **R. LA BRIE**

JOB NO. **11-1018**

DATE **4/7/10**


SHEET

**4**

OF **4** SHEETS

SEISMIC ANCHORAGE

ELEVATED FLOOR



**NOTES:**

1. ANCHORAGE DESIGN PER 2010 CALIFORNIA BUILDING CODE - SECTION 1613A AND ASCE 7-05 SECTIONS 12 AND 13. ALLOWABLE STRESS DESIGN IS USED.

HORIZONTAL FORCE ( $E_h$ ) =  $0.97 W_p$  ( $S_{Ds} = 1.93$ ,  $a_p = 1.0$ ,  $I_p = 1.5$ ,  $R_p = 2.5$ )

VERTICAL FORCE ( $E_v$ ) =  $0.27 W_p$

$T_{MAX} = 262$  LBS/BOLT  
 $V_{MAX} = 240$  LBS/BOLT

2. CENTER OF GRAVITY (C.G.) WEIGHT IS A MAXIMUM. THIS PRE-APPROVAL ENCOMPASSES ALL WEIGHTS UP TO THE MAXIMUM WEIGHT SHOWN.
3. ENGINEER OF RECORD FOR THE BUILDING SHALL PROVIDE SUPPORT STRUCTURE DESIGNED TO SUPPORT WEIGHTS AND FORCES SHOWN, IN ADDITION TO ALL OTHER LOADS.
4. SEE GENERAL NOTES: SHEETS 1-2

