



A

TRANSPORT CART  
PLAN VIEW  
1/2" = 1"

B

TRANSPORT CART  
ISOMETRIC VIEW  
1/2" = 1"

CART WEIGHT	83.75 lbs
CART MAX LOAD	379.25 lbs



C

TRANSPORT CART  
FRONT VIEW  
1/2" = 1"

D


TRANSPORT CART  
SIDE VIEW  
1/2" = 1"

## CA969 CART

CART FOR CA9618-HS1 / HS2 STERILIZER


TECHNICAL DRAWING

Dwg. No. :	CUT-CA969 CART	Rev. :	A
Drawn By :	A. HERRON	Date :	6/28/16
Scale :	1/2" = 1"	Page :	1 of 2




**LOADING RACK  
PLAN VIEW**  
1/2" = 1"

A



**LOADING RACK  
FRONT VIEW**  
1/2" = 1"

B



**LOADING RACK  
SIDE VIEW**  
1/2" = 1"

C

RACK WEIGHT	39.5 lbs
RACK MAX LOAD	299.75 lbs
SHELF WEIGHT	13.25 lbs
SHELF MAX LOAD	99.25 lbs

## CA969 RACK

RACK FOR CA9618-HS1 / HS2 STERILIZER

TECHNICAL DRAWING

Dwg. No. :	CUT-CA969 CART	Rev. : A
Drawn By :	A. HERRON	Date : 6/28/16
Scale :	1/2" = 1'	Page : 2 of 2