

# Enzymatic Pretreatment Foam Spray



## Specially formulated multi-enzyme dual foam/liquid pretreatment spray

For pre-cleaning and transportation of surgical grade instruments, reusable items, scopes and other devices

The Belimed Protect® Enzymatic Pretreatment Foam Spray is designed to offer superior cleaning performance at point of use. The triple enzyme formula initiates cleaning upon contact, providing increased effectiveness in removing blood, fibrin, fats and other bioburden. The pretreatment can be used as a foam or liquid spray and is designed to be compatible with Belimed Protect solutions for manual cleaning, ultrasonic and automated wash processes. Belimed Protect products also have a unique color and number coded labeling system designed to reinforce what product is used at each step and reduce the risk of dosing errors.

Subtilisin is a common protease that has been associated with allergies, including skin and eye irritation and respiratory reactions. The Belimed Protect enzymatic products are subtilisin-free, providing a better work environment for CSSD staff. Additionally, all products have been independently tested for dermal and ocular irritation, and confirmed to be non-irritating to skin and eyes.

### Product Features

- pH Neutral
- Gentle on instruments
- Easily rinsed
- Fragrance-free
- Stays wet up to 72 hours
- Safe for end users
- Reduces risk of biofilm formation

### Key Formula Ingredients

- Three (3) enzymes (protease for proteins, lipase for fatty acids and amylase for starches)
- Unique blend of surfactants solubilize bioburden and retain wetness
- Corrosion inhibitor

Key Chemical Attributes	
pH	6.5-8.5
Density	1.010-1.020
Color	Blue

### Material Compatibility

The Belimed Protect Enzymatic Pretreatment Foam Spray demonstrates a wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- Aluminum
- Anodized Aluminum
- Titanium
- Copper
- Brass
- All commonly used plastics and elastomers, including Polyurethane, Silicone and Teflon
- Synthetics (Vinyl and Nitrile)

### Storage

Do not store product at extreme hot or cold temperatures

### Safety Information

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

### Recommended Practice:

- After each case, staff should spray all instrument trays
- The instruments should be covered completely with product
- After a few hours, check trays and re-spray if necessary

### Instructions for Use

1. Shake bottle before use
2. Spray instruments and equipment as soon as possible after procedure
3. Keep instruments and equipment wetted with solution until decontamination process is started
4. Follow facility policies and procedures

### Product Ordering Information

	Color Code	Item Order #	Product ID (On the Label)	Pkg	Units per case
Enzymatic Pretreatment Foam Spray	1	106-8088	2110032	32 oz	12
		1008899	2111002	2 oz	48
		110-2632	2110001	1 oz	4
		110-2633	2110002	2.5 oz	2