

#### NOTES:

##### 1. MATERIAL:

SINK BASIN: 14 GAUGE  
304 STAINLESS STEEL  
SINK SKIRT: 16 GAUGE  
304 STAINLESS STEEL


2. ALL PLUMBING 1/2 COPPER OR BRASS PIPE.  
THE SINK WILL BE FURNISHED WITH IN LINE CHECK VALVES ON SUPPLY LINES.

3. UTILITY REQUIREMENT  
COLD WATER 20 TO 50 PSIG 70 DEG F MAX  
HOT WATER 20 TO 50 PSIG 120 DEG F TO 140 DEG F

4. POWER REQUIREMENTS (FOR SENSOR OPERATIONS)  
120 VAC, 60 HZ, SINGLE PHASE, 3.0 AMP OR  
220 VAC, 60 HZ, SINGLE PHASE, 1.5 AMP.  
INSTALLED PER LOCAL BUILDING CODES.

5. NET WEIGHT OF SINK 130 LB..

## SINGLE SCRUB SINK SS32



820 South Mulberry  
Millsstadt, IL 62260  
877-828-9975

NOTES:


1. MATERIAL:

SINK BASIN: 14 GAUGE  
304 STAINLESS STEEL  
SINK SKIRT: 16 GAUGE  
304 STAINLESS STEEL

2. ALL PLUMBING 1/2" COPPER OR BRASS PIPE.  
THE SINK WILL BE FURNISHED WITH IN LINE CHECK  
VALVES ON SUPPLY LINES.

3. UTILITY REQUIREMENT  
COLD WATER 20 TO 50 PSIG 70 DEG F MAX  
HOT WATER 20 TO 50 PSIG 120 DEG F TO 140 DEG F

4. POWER REQUIREMENTS (FOR SENSOR OPERATIONS)  
120 VAC, 60 HZ, SINGLE PHASE, 3.0 AMP OR  
220 VAC, 60 HZ, SINGLE PHASE, 1.5 AMP.  
INSTALLED PER LOCAL BUILDING CODES.




# SINGLE SCRUB SINK PEDESTAL BASE

SS32-PB




820 South Mulberry  
Millsstadt, IL 62260  
877-828-9975



SINGLE SCRUB SINK  
WALL MOUNTING  
BRACKETS INCLUDED WITH SINK  
SS32

**mac**  
MEDICAL

820 South Muberry  
Millstadt, IL 62260  
877-828-9975



SINGLE SCRUB SINK  
RECOMMENDED OPTIONAL IN WALL  
CHAIR CARRIER SYSTEM  
SS32



820 South Mulberry  
Millstadt, IL 62260  
877-828-9975