



## Belimed Protect® Corner Guards



### Effective, Efficient and Environmentally-Friendly

A unique solution to prevent rips and tears in wrapped surgical trays

#### Effective

Belimed Protect® Corner Guards are validated for use with steam, ethylene oxide and plasma sterilization. The patented material will not retain moisture, unlike towel or foam alternatives. Specially designed for large surgical and loaned trays, they will not flatten under heavy weight and protect wrap from puncturing tray feet. The extended edges ensure that the corner guards stay in place through sterilization, storage and transport.

#### Environmentally-Friendly

Our unique resin material is latex-free and non-toxic. Belimed Protect® Corner Guards can be reused multiple times and are less wasteful than single-use solutions. They are also recyclable to further reduce environmental impact.


#### Efficient

Belimed Protect® Corner Guards prevent tears, holes and rips in sterilization wrap by creating a barrier from tray feet. Eliminating holes saves time, money and frustration associated with rework. The guards are reusable, meaning they can be used up to 10 times, making them the most cost effective way to prevent tears in sterile wrap.





## Features.



## Material Specifications.

- Thermoplastic elastomer
- Latex-Free
- Patented
- FDA Compliant
- Made in USA



## Ordering Information

Product Specifications	Order#	Qty
Belimed Protect® Corner Guard, Case	111-2840	500
Belimed Protect® Corner Guard, Sleeve	111-2839	100