



Belimed

Belimed | Protect®

Lubricant WD

Enzymatic Detergent WD

5

Instrument Lubricant for use in
Automated Washer Disinfectors

4

Multi-Enzyme, subtilisin-free neutral detergent formulated for
Automated Washer Disinfectors
Concentrate, No Foam

2.5 Gallons (9.4L)

2.5 Gallons (9.4L)

Belimed | Protect®



ABSOLUTE Confidence



WE ARE DIFFERENT!

Our singular focus on your world gives us the expertise to help you uncover and resolve process challenges, streamline and improve the overall work environment and provide the OR with safe and sterile instruments.

We combine our deep knowledge of sterile processing with unique product and service solutions to drive out cost, improve efficiency, and deliver quality.

That's Sterile Processing with Absolute Confidence.

Belimed Protect®

Proprietary, exclusive 3-enzyme formulation designed and developed by Belimed:

- Optimizes the cleaning process delivered by the superior low impingement mechanical action of Belimed's washer/disinfectors combined with our specially formulated cleaning solutions.
- Minimizes the variability associated with competitive detergents not designed for Belimed equipment - which can lead to higher cost and cycle times.
- Reduces dosing requirements for manual and ultrasonic cleaning.
- Ensures solution compatibility throughout the cleaning process, enhancing cleaning efficacy.

Subtilisin-free formulation:

- May reduce risk of allergies and respiratory issues, creating a safer work environment for CSSD staff.
- May contribute to decreasing the risk of TASS (Toxic Anterior Segment Syndrome) associated with enzyme use in ophthalmic instrument reprocessing.

Color coded packaging designed to correspond to proper equipment:

- Easy to locate and connect detergents and equipment quickly and efficiently.
- May reduce risk from:
 - Foaming and faulting of equipment
 - Potential rework due to incorrect dosing
 - Faulty wash monitor results

At A Glance

- **Multi-Enzyme Formulation Specially Designed For Each Cleaning Step**
- **Proprietary Subtilisin-Free Formulation**
- **Fragrance-Free**
- **pH Neutral**
- **Biodegradable**
- **Gentle/Non-Corrosive**
- **Wide Range of Material Compatibility**
- **Color Coded**



Belimed Protect® Cleaning Solutions

Belimed's Protect line of proprietary cleaning solutions reinforces our commitment to high performance Sterile Processing Environments. We focus on what's most important - Protecting your patients, staff, environment and your bottom line. The Protect portfolio is specially formulated for all steps of the decontamination process.



Multi-Enzyme Instrument Pre-Treatment Foam Spray



Multi-Enzyme Detergent for Manual Cleaning



Multi-Enzyme Detergent for Ultrasonic Washers

Solutions for Pre-Treatment, Manual and Ultrasonic Cleaning

Exclusive multi-enzyme formulations designed by Belimed for instrument pre-treatment, manual and ultrasonic cleaning as part of our family of solutions for all cleaning processes. The subtilisin free formulation reduces odors and risk of allergic reactions, creating a better environment for staff.

* *Protease, lipase and amylase*



Multi-Enzyme Detergent for Automated Washer/Disinfectors (Concentrate)



Multi-Enzyme Detergent for Automated Washer/Disinfectors (Concentrate-PLUS)



Lubricant for Automated Washer/Disinfectors

Solutions for Automated Washer/Disinfectors

Proprietary multi-enzyme formulations designed by Belimed to optimize our low impingement wash process, improving equipment cleaning performance. Our enzyme and lubricant solutions deliver consistently clean instruments. The subtilisin free formulation reduces odors and risk of allergic reactions, creating a better environment for staff. Products demonstrate equal efficacy in high impingement washers.



Neutral Non-Enzymatic Detergent for Cart Washers and Automated Washer/Disinfectors



Cart Rinse

Solutions for Cart Washers

Specially formulated solutions designed for the automated cart cleaning process. Our neutral non-enzymatic and cart rinse products assist in quick and effective cleaning and turnaround of surgical case carts in the busy CSSD. Concentrated formula reduces need for frequent solution replacement. Color coded labels make it easy for cart washer setup.



Ready-to-Use Instrument Lubricant



Descaler/Delimer



Rust and Stain Remover



Adhesive Remover



Alkaline Detergent



Neutralizer Detergent

Ancillary Items

Ancillary items complete our portfolio, providing you additional products that support efficient instrument preparation and preservation, and maintenance of instruments, washers and sterilizers.

Accessories



Automated Dosing Pump



Proportioner



Wall Mount Shelf

Product Descriptions

Color Code	Numeric Code**	Item #	Product Description	Size	Product Use
	1	2110032	Multi Enzyme Instrument Pre-Treatment Foam Spray	32 oz, 1 gal, 2.5 gal	Instrument Pre-Treatment and Transport
	2	21200XX	Multi Enzyme Detergent for Manual Soaking and Cleaning, Concentrate, Low Foam	1 gal, 2.5 gal, 5 gal, 15 gal	Manual Washing/Sink
	2	21210XX	Multi Enzyme Detergent for Manual Soaking and Cleaning, Concentrate, Low Foam, with Fragrance	1 gal, 2.5 gal, 5 gal	Manual Washing/Sink
	3	21300XX	Multi Enzyme Detergent for Ultrasonic Cleaning, Ultra Concentrate, No Foam	1 gal, 2.5 gal, 5 gal	Ultrasonic Washers
	3	21310XX	Multi Enzyme Detergent for Ultrasonic Cleaning, Ultra Concentrate, No Foam, with Fragrance	1 gal, 2.5 gal, 5 gal	Ultrasonic Washers
	4	21400XX	Multi Enzyme Detergent for Automated Washer/Disinfectors, Concentrate, No Foam	1 gal, 2.5 gal, 5 gal, 15 gal, 30 gal	Automated Washer/Disinfectors
	4+	21410XX	Multi Enzyme Detergent for Automated Washer/Disinfectors, Concentrate PLUS, No Foam	1 gal, 2.5 gal, 5 gal	Automated Washer/Disinfectors
	5	21500XX	Lubricant for Automated Washer/Disinfectors	1 gal, 2.5 gal, 5 gal, 15 gal, 30 gal	Automated Washer/Disinfectors
	6	21600XX	Neutral Non-Enzymatic Detergent for Cart Washers and Automated Washer/ Disinfectors, Concentrate, Low Foam	1 gal, 2.5 gal, 5 gal, 15 gal, 30 gal, 55 gal	Automated Cart Washers, Automatic Washer/Disinfectors, Manual Washing Sink
	7	21700XX	Cart Rinse	1 gal, 2.5 gal, 5 gal, 15 gal, 30 gal, 55 gal	Automated Cart Washers
		2150032	Ready to Use Instrument Care Lubricant	32 oz	Prep and Pack
		2181034	Descaler/Delimer	1 gal, 5 gal	Automated Cart Washers, Automatic Washer/Disinfectors
		2182034	Adhesive Remover	8 oz	Decontamination
		2183034	Rust and Stain Remover	1 gal	Decontamination
		219100X	Alkaline Detergent	2.5 gal, 5 gal	Automated Washer/Disinfectors
		219200X	Neutralizing Detergent	2.5 gal, 5 gal	Automated Washer/Disinfectors, Decontamination
		BMP-100	Manual Dose Pump for Sink	Each	
		BPP-100	Proportioner for Sink	Each	
		BDP-100	Automated Dose Pump for Sinks	Each	
		BWM-100	Wall Mount for Detergent Storage	Each	

**Numeric code provided as guide for color blindness



Belimed
Infection Control

1535 Hobby Street, Suite 103
North Charleston, SC 29405

1.800.451.4118
1.843.216.7424

info.us@belimed.com

201601_Protect_Rev. 02

