

Cart Rinse



Rinse Aid for Cart Washers and Automated Washer/Disinfectors

Rinse Aid for spot-free case carts and other stainless steel products.

Product Features:

- Ensures complete drying of case carts in automated cart washers and other automated equipment
- Conditions water for spot free drying and sheeting action reduces dry time
- Gentle ingredients will not damage case carts, casters or containers
- Biodegradable and aquatic life safe
- pH neutral and non-corrosive
- Fragrance free

Key Chemical Attributes

pH Neat	1.5-3.5
pH Use Dilution	5.0-7.0
Density	1.02-1.03
Color	Dark Blue

Handling and Storage:

Do not store at extreme hot or cold temperatures.


Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

Instructions for Use:

1. Add 1/8 to 1/2 oz of rinse aid (1 to 4 ml/L) per gallon of water.
2. Follow facility policies and procedures.

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Cart Rinse		106-8079	2170002	2.5 gal	2
		106-8080	2170005	5 gal	1
		106-8081	2170015	15 gal	1
		106-8082	2170030	30 gal	1
		106-8083	2170055	55 gal	1



For more information, visit belimed.com

REV12.22.2020_Q0

Engineers of Confidence.