



Processing Sink Pegboard and Accessories



Pegboards and accessory options, for both fixed and adjustable height sinks, optimize the functionality of your sink and improves efficiency and workflow.



Processing Sink Pegboard and Accessories



- 1 Pegboard** - ¼" holes on 1" centers
 - 120" W x 36" H for Triple Bay Sinks
 - 96" W x 36" H for Double Bay Sinks
 - 72" W x 36" H for Single Bay Sinks
- 2 LED Overhead Light**
- 3 Glovebox Holder** - Holds 2 boxes
 - 10" W x 3" D x 10" H
- 4 Bottle Holder** - Holds 2 bottles, up to 3 5/8" diameter
 - 9" W x 5" D x 5" H
- 5 Gallon Container Holder** - Holds 2 containers
 - 13" W x 6 1/2" D x 5" H
- 6 Brush Holder** - Includes 10 hooks
 - 9 1/4" W x 1" D x 1 1/4" H
- 7 Small Bin**
 - 4 1/8" W x 7 3/8" D x 3" H
- 8 Large Bin**
 - 5 1/2" W x 10 7/8" D x 5" H
- 9 Display Arm and Mount**
 - Holds up to a 25 lb Display
 - Display not included
- 10 Shelf**
 - 15" D x Various Lengths
- 11 Bin Rail**
 - Used for holding small and large bins - 29" W

REV4.5.2021_Q320

Engineers of Confidence.