




WA 290 Automation System

Automated Loading and unloading for the
WD 290 IQ washer-disinfector

Engineers of Confidence.

Belimed
Infection Control



WA 290 Automation System

Designed to make the Decontamination side more effective and efficient with loading and washing.


Simple, efficient, automated - Load the rack and the WA 290 does the rest

- Just load the rack on the feed line
- Shuttle will pick up a rack
- Deliver it to the nearest open washer
- Shuttle includes sensors to avoid collisions
- Temporarily suspend the shuttle activity to manually load a washer

Easy to use and configurable for your needs

- Single-point loading
- Various rack queue quantities and configurations to meet processing needs and available space
- Keep the racks on the system with the Cross Conveyor or
- Remove racks using a cart and Direct Dock

Feed Line Queuing Options



Designed to minimize the footprint on the clean side and provide process flexibility.

Optimizes floor space and improves visibility to the finished racks

- Requires very little of your valuable clean side space
- Racks travel along the washers on the conveyor system
- Once the racks are unloaded send them back through the automated gate using the return line

Easy to use and adaptable to your process

- The transverse conveyor system provides a clear view of all racks
- Keep the racks on the system and unload the trays to a cart or
- Remove the entire rack, onto a cart, once it travels to the unload station

Contact

Belimed Inc. is a leading supplier of medical and surgical instrument sterilization, disinfection, and cleaning products and services.

With 12 Belimed companies throughout Europe, North America and China, Belimed is represented worldwide by a strong network in over 80 countries.

USA Headquarters

Belimed, Inc.
8351 Palmetto Commerce Parkway
Suite 101
Ladson, SC 29456
1.800.451.4118
info.us@belimed.com

Global Headquarters

Belimed AG
Grienbachstrasse 11
6300 Zug
Switzerland
+41 41 449 78 88
info@belimed.com

REV2.7.2022_Q320



Engineers of Confidence.

Belimed
Infection Control

For more information,
visit www.belimed.com

